

PENPREP EVO™

Cartridge-based reconstitution


PENPREP **EVO** incorporates patented vacuum transfer technology, which enables an accurate reconstitution process, while minimizing both drug hold-up in the vial and air transfer back to the cartridge.

The modular concept of PENPREP EVO ensures it can easily be configured for 13mm or 20mm drug vials and optimized for different diluent fill volumes.

Simplify your Pen injector cartridge preparation with PENPREP **EVO™**.



DUOJECT
MEDICAL
SYSTEMS



DUOJECT.COM