

# INTER-VIAL™

## Pre-assembled cartridge-based reconstitution

The drug is contained in a standard drug vial and the diluent in a standard ISO cartridge. Both containers are pre-installed inside the device minimizing the steps for reconstitution and drug delivery.

Whether you need parenteral or topical application, our patented activation system allows an intuitive transfer to reconstitute drug and withdraw admixture to the detachable delivery unit.

Differentiate your product with INTER-VIAL™



**DUOJECT**  
MEDICAL  
SYSTEMS

50 chemin de Gaspé, Complex B-5, Bromont (Québec), Canada J2L 2N8 T: (450) 534-3666



DUOJECT.COM